

Product datasheet

Recombinant Human Density Regulated Protein ab111640

1 图像

概述

产品名称	重组人Density Regulated蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

Accession	<a href="#">O43583</a>
种属	Human
序列	<b>MGSSHHHHHSSGLVPRGSHMAADISESSGADCKGDPRNSAKLDADYPLR VLYCGVCSLPTEYCEYMPDVAKCRQWLEKNFPNEFAKLTVENSPKQEAGI SEGQGTAGEEEEEKKQKRGGRGQIKQKKKTVPQKVTIAKIPRAKKKYVTR VCGLATFEIDLKEAQRFFAQKFSCGASVTGEDEIIIQGDFTDDIIDVIQE KWPEVDDDSIEDLGEVKK</b>
分子量	24 kDa including tags
氨基酸	1 to 198
标签	His tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab111640** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Mass Spectrometry
质谱法	MALDI-TOF
纯度	> 85 % SDS-PAGE. ab111640 is purified using conventional chromatography.
形式	Liquid

## 稳定性和存储

Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.

pH: 8.00

Constituents: 20% Glycerol, 0.04% DTT, 0.88% Sodium chloride

Note: Contains 89% Tris-HCl buffer.

## 常规信息

### 功能

May be involved in the translation of target mRNAs by scanning and recognition of the initiation codon. Plays a role in the modulation of the translational profile of a subset of cancer-related mRNAs when recruited to the translational initiation complex by the oncogene MCTS1.

### 组织特异性


Highly expressed in heart and skeletal muscle and moderately expressed in the brain, placenta, liver and pancreas. Weakly expressed in the lung and kidney.

### 序列相似性

Belongs to the DENR family.

Contains 1 SU1 domain.

## 图片



15% SDS-PAGE analysis of ab111640 (3 µg)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors