

## Product datasheet

# Recombinant Human Galectin 10 protein ab117473

### 概述

---

|      |                  |
|------|------------------|
| 产品名称 | 重组人Galectin 10蛋白 |
| 蛋白长度 | Protein fragment |

### 描述

---

|    |                  |
|----|------------------|
| 性质 | Recombinant      |
| 来源 | Escherichia coli |

### 氨基酸序列

Accession [Q05315](#)

种属 Human

序列  
 MGSSHHHHHH SGLVPRGSH MSLLPVPYTE  
 AASLSTGSTV TIKGRPLACF LNEPYLQVDF  
 HTEMKEESDI VFHFQVCFGR RVVMNSREYG  
 AWKQQVESKN MPFQDQEFEL SISVLPDKY  
 QVMVNGQSSY TFDHRIKPEA VKMVQVWRDI  
 SLTKFNVSYL KR

|     |                       |
|-----|-----------------------|
| 分子量 | 19 kDa including tags |
| 氨基酸 | 1 to 142              |
| 标签  | His tag N-Terminus    |

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab117473** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

|    |        |
|----|--------|
| 形式 | Liquid |
|----|--------|

### 制备和贮存

---

|        |   |
|--------|---|
| 稳定性和存储 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.<br>pH: 7.50<br>Constituents: Tris acetate, 20% Glycerol, 0.03% EDTA |
|--------|---|

## 常规信息

---

|       |   |
|-------|---|
| 功能    | May have both lysophospholipase and carbohydrate-binding activities.  |
| 组织特异性 | Expressed exclusively by eosinophils and basophils. Not detected in monocytes and neutrophils.                          |
| 序列相似性 | Contains 1 galectin domain.   |
| 细胞定位  | Cytoplasmic granule. Localized in granules from where it may be secreted or transported to other locations in the cell. |

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors