

Recombinant Human HSC70 Interacting Protein HIP ab61296

1 图像

描述

产品名称	重组人HSC70 Interacting蛋白HIP
纯度	> 95 % SDS-PAGE. >99% by SDS-PAGE
表达系统	Escherichia coli
蛋白长度	Full length protein
无动物成分	No
性质	Recombinant
种属	Human
描述	重组人HSC70 Interacting蛋白HIP

技术指标

Our **Abpromise guarantee** covers the use of **ab61296** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Western blot
形式	Liquid
补充说明	This product is manufactured by BioVision, an Abcam company and was previously called 4904 HIP, human recombinant. 4904-25 is the same size as the 25 µg size of ab61296.

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot. Store at -80°C. Avoid freeze / thaw cycle. pH: 7.20 Constituents: 0.418% MOPS, 1.875% Potassium chloride, 0.0078% 2-Mercaptoethanol, 0.0146% EDTA
--------	--

常规信息

功能	One HIP oligomer binds the ATPase domains of at least two HSC70 molecules dependent on
----	--

activation of the HSC70 ATPase by HSP40. Stabilizes the ADP state of HSC70 that has a high affinity for substrate protein. Through its own chaperone activity, it may contribute to the interaction of HSC70 with various target proteins.


序列相似性

Belongs to the FAM10 family.
Contains 1 ST11 domain.
Contains 3 TPR repeats.

细胞定位

Cytoplasm.

图片



Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors