

## Product datasheet

# Recombinant Human IFIT2 protein ab126008

### 概述

---

产品名称	Recombinant人IFIT2 protein
蛋白长度	Protein fragment

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession [P09913](#)

种属	Human
分子量	26 kDa
氨基酸	100 to 321

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab126008** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilised

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Store at -20°C. Constituents: 0.32% Tris HCl, 0.58% Sodium chloride
复溶	Reconstitute with water to desired concentration.

### 常规信息

---

序列相似性	Belongs to the IFIT family. Contains 6 TPR repeats.
-------	--

---

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors