

Product datasheet

Recombinant Human KIAA0101 protein ab126015

概述

产品名称	Recombinant人KIAA0101 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	Q15004
种属	Human
分子量	10 kDa
氨基酸	8 to 101

技术指标

Our [Abpromise guarantee](#) covers the use of **ab126015** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilised
补充说明	50ug package, it is suggested to dissolve directly in 50ul DDW and the final concentration is 1.0ug/ul

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Constituents: 0.31% Tris HCl, 0.58% Sodium chloride
复溶	Reconstitute with water to desired concentration.

常规信息

功能	May be involved in protection of cells from UV-induced cell death.
组织特异性	Expressed predominantly in liver, pancreas and placenta. Not detected in heart or brain. Highly

expressed in a number of tumors, especially esophageal tumors, in anaplastic thyroid carcinomas and in non-small-cell lung cancer lines.

细胞定位

Nucleus. Mitochondrion.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors