

## Product datasheet

# Recombinant Human LECT1 protein ab112725

### 概述

---

产品名称	重组人LECT1蛋白
蛋白长度	Protein fragment

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession [O75829](#)

种属	Human
分子量	38 kDa

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab112725** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
形式	Lyophilised

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. pH: 6.80 Constituents: 1.58% Tris HCl, 3.09% DTT, 4% SDS, 0.2% Bromophenol blue, 20% Sucrose
复溶	Add 100ul of distilled water.

### 常规信息

---

功能	Bifunctional growth regulator that stimulates the growth of cultured chondrocytes in the presence of basic fibroblast growth factor (FGF) but inhibits the growth of cultured vascular endothelial cells. May contribute to the rapid growth of cartilage and vascular invasion prior to the
----	--

replacement of cartilage by bone during endochondral bone development. Inhibits in vitro tube formation and mobilization of endothelial cells. Plays a role as antiangiogenic factor in cardiac valves to suppress neovascularization.

**组织特异性**

Detected in cartilage and cardiac valves (at protein level). Detected in the laminae fibrosa, spongiosa and ventricularis layers of normal cardiac valves (at protein level). Expression is decreased cardiac valves of patients with valvular heart disease (at protein level). Weakly expressed in chondrosarcoma.

**序列相似性**

Belongs to the chondromodulin-1 family.  
Contains 1 BRICHOS domain.

**发展阶段**

Expressed at 9 weeks in developing cartilagenous bone rudiments.

**翻译后修饰**

After cleavage, the post-translationally modified ChM-I is secreted as a glycoprotein.

**细胞定位**

Endomembrane system and Secreted. Accumulated in the inter-territorial matrix of cartilage.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors