

## Product datasheet

# Recombinant Human Liver Arginase protein ab117182

### 概述

产品名称	重组人Liver Arginase蛋白
蛋白长度	Full length protein

### 描述

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession [P05089](#)

种属 Human

序列  
 MSAKSRITIGI IGAPFSKGQP RGGVEEGPTV  
 LRKAGLLEKL KEQECDVKDY GDLPFADIPN  
 DSPFQIVKNP RSVGKASEQL AGKVAEVKKN  
 GRISLVLGGD HSLAIGSISG HARVHPDLGV  
 IWVDAHTDIN TPLTTTSGNL HGQPVSFLLK  
 ELKGKIPDVP GFSWVTPCIS AKDIVYIGLR  
 DVDPGEHYIL KTLGIKYFSM TEVDRLGIGK  
 VMEETLSYLL GRKKRPIHLS FDVDGLDPSF  
 TPATGTPVVG GLTYREGLYI TEEIYKTGLL  
 SGLDIMEVNP SLGKTPEEVT RTVNTAVAIT  
 LACFGLAREG NHKPIDYLNPKLEHHHHHHH

分子量 36 kDa including tags

氨基酸 1 to 322

标签 His tag C-Terminus

额外的序列信息 ARG1 Human Recombinant produced in E.Coli is a single, non-glycosylated, polypeptide chain containing 330 amino acids (1-322a.a.) and having a molecular mass of 35.8kDa.

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab117182** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用 SDS-PAGE

<b>纯度</b>	> 85 % SDS-PAGE. ab117182 is purified by standard chromatography.
<b>形式</b>	Liquid
<b>补充说明</b>	Although stable at 4°C for 1 week, ab117182 should be stored desiccated below -18°C. Please prevent freeze thaw cycles.

## 制备和贮存

---

<b>稳定性和存储</b>	Shipped at 4°C. Please see notes section. pH: 8.00 Constituents: 0.32% Tris HCl, 0.03% DTT, 0.58% Sodium phosphate, 20% Glycerol
---------------	--

## 常规信息

---

<b>通路</b>	Nitrogen metabolism; urea cycle; L-ornithine and urea from L-arginine: step 1/1.
<b>疾病相关</b>	Defects in ARG1 are the cause of argininemia (ARGIN) [MIM:207800]; also known as hyperargininemia. Argininemia is a rare autosomal recessive disorder of the urea cycle. Arginine is elevated in the blood and cerebrospinal fluid, and periodic hyperammonemia occurs. Clinical manifestations include developmental delay, seizures, mental retardation, hypotonia, ataxia, progressive spastic quadriplegia.
<b>序列相似性</b>	Belongs to the arginase family.
<b>细胞定位</b>	Cytoplasm.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors