

Product datasheet

Recombinant Human PMP2 protein ab83017

1 图像

概述

产品名称	重组人PMP2蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
种属	Human

技术指标

Our [Abpromise guarantee](#) covers the use of **ab83017** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 95 % SDS-PAGE.
形式	Liquid


制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C. Preservative: None Constituents: 50% Glycerol, PBS, pH 7.2
--------	---------------------------------------------------------------------------------------------------------------------------------

常规信息

功能	May play a role in lipid transport protein in Schwann cells. May bind cholesterol.
序列相似性	Belongs to the calycin superfamily. Fatty-acid binding protein (FABP) family.
结构域	Forms a beta-barrel structure that accommodates hydrophobic ligands in its interior.

## 图片



SDS-PAGE showing ab83017.

SDS-PAGE - PMP2 protein (His tag) (ab83017)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors