

Product datasheet

Recombinant Human PTP mu protein ab71346

1 图像

概述

产品名称	Recombinant人PTP mu protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Baculovirus infected Sf9 cells
氨基酸序列	
种属	Human
氨基酸	764 to 1452

技术指标

Our [Abpromise guarantee](#) covers the use of **ab71346** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Western blot
纯度	> 90 % SDS-PAGE. Purity was determined to be >90% by densitometry.
形式	Liquid

制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 25% Glycerol, 50mM Tris HCl, 150mM Sodium chloride, 0.25mM DTT, 0.1mM PMSF, pH 7.5
--------	---

常规信息

功能	Involved in cell-cell adhesion through homophilic interactions. May play a key role in signal transduction and growth control.
序列相似性	Belongs to the protein-tyrosine phosphatase family. Receptor class 2B subfamily. Contains 4 fibronectin type-III domains.

Contains 1 Ig-like C2-type (immunoglobulin-like) domain.
Contains 1 MAM domain.
Contains 2 tyrosine-protein phosphatase domains.

细胞定位

Membrane.

Recombinant Human PTP mu protein 图像



SDS-PAGE showing ab71346 at approximately 110kDa.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors