

Product datasheet

Recombinant Human PURA protein ab114480

1 图像

概述

产品名称	Recombinant人PURA protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ
氨基酸序列	
Accession	Q00577
种属	Human
序列	TQGQTIALPAQGLIEFRDALAKLIDDYGVVEEPAELPEGTSLTVDNKRFF FDVGSNKYGVFMRVSEVKPTYRNSITVPYKVVAKFGHTFCKYSEEMKKIQ EKQREKRAAC
分子量	38 kDa including tags
氨基酸	183 to 292

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114480** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml. Best used within three months from the date of receipt.


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00 Constituents: 0.3% Glutathione, 0.79% Tris HCl
--------	---

常规信息

功能	This is a probable transcription activator that specifically binds the purine-rich single strand of the PUR element located upstream of the MYC gene. May play a role in the initiation of DNA replication and in recombination.
序列相似性	Belongs to the PUR DNA-binding protein family.
细胞定位	Nucleus.

Recombinant Human PURA protein 图像



12.5% SDS-PAGE showing ab114480 at approximately 37.73kDa stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors