

Product datasheet

Recombinant Human ROM1 protein ab116995

1 图像

概述

产品名称	Recombinant人ROM1 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession [Q03395](#)

种属 Human

序列 QLRyhccGRHGyKdWFGVQwVSSRYLDPGDRDVADRIQSNVEGLYLTdGV
PFSCCNPHSPRPCLQNRLSDSYAHPLFDPRQPNQNLWAQGCHEVLLLEHLQ
DL

分子量 37 kDa including tags

氨基酸 163 to 264

技术指标

Our [Abpromise guarantee](#) covers the use of **ab116995** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	SDS-PAGE
	Western blot

形式 Liquid

补充说明 Protein concentration is above or equal to 0.05 mg/ml.
This protein is best used within three months from the date of receipt.

制备和贮存


稳定性和存储 Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00
Constituents: 0.79% Tris HCl, 0.3% Glutathione

常规信息

功能	May function as an adhesion molecule involved in stabilization and compaction of outer segment disks or in the maintenance of the curvature of the rim. It is essential for disk morphogenesis.
组织特异性	Retina (photoreceptor). In rim region of ROS (rod outer segment) disks.
疾病相关	Defects in ROM1 may cause retinitis pigmentosa (RP) [MIM:268000]; when associated with defects in PRPH2.
序列相似性	Belongs to the PRPH2/ROM1 family.
细胞定位	Membrane.

图片



12.5% SDS-PAGE showing ab116995 at
approx. 36.85 kDa.
Stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors