

Product datasheet

Recombinant Human SNRPC protein ab116159

1 图像

概述

产品名称	Recombinant人SNRPC protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	P09234
种属	Human
序列	MGSSHHHHHHSSGLVPRGSHMGSMKPKFYCDYCDTYLTHDSPSVRKTHCSG RKHKENVKDYYQKWMEEQAQSLIDKTTAAAFQQGKIPPTPFSAPPPAGAMI PPPSLPGPPRPGMMPAPHMGGPPMMPMMGPPPPGMPVGPAP GMRPP MGGHMPMPGPPMMRPPARPMMPVTRPGMTRPDR
分子量	20 kDa including tags
氨基酸	1 to 159
标签	His tag N-Terminus

技术指标

Our [Abpromise guarantee](#) covers the use of **ab116159** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
质谱法	MALDI-TOF
纯度	> 85 % SDS-PAGE. purified by using conventional chromatography
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
--------	--


pH: 8.00

Constituents: 0.32% Tris HCl, 50% Glycerol, 1.75% Sodium chloride, 0.08% DTT, 5.84% EDTA

常规信息

功能	This protein is associated with snRNP U1.
序列相似性	Belongs to the U1 small nuclear ribonucleoprotein C family. Contains 1 matrin-type zinc finger.
细胞定位	Nucleus.

Recombinant Human SNRPC protein 图像



15% SDS-PAGE showing ab116159 at approximately 19.8kDa (3 μ g)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors