

## Product datasheet

# Recombinant Human SS18 protein ab116825

### 1 图像

#### 概述

---

产品名称	Recombinant人SS18 protein
蛋白长度	Full length protein

#### 描述

---

性质	Recombinant
来源	Wheat germ

#### 氨基酸序列

Accession [Q15532](#)

种属 Human

序列 MSVAFAAPRQRGKGEITPAAIQKMLDDNNHLIQCIMDSQNKGKTSECSQY  
 QQMLHTNLVYLATIADSNQNMQSLLPAPPTQNMMPMGPGGMNQSGLPPPPR  
 SHNMPSDGMVGGGPPAPHMQNMNGQMPGPNHMPMQGPGPNQLNMTNSSM  
 NMPSSSHGSMGGYNHVPSSQSMPVQNQMTMSQGGQPMGNYGPRPNMSMQP  
 NQGPMMHQQPPSQYNMPQGGGQHYQGQQPPMGMMQVQNGNHMMGQRQI  
 PPYRPPQQGPPQQYSGQEDYYGDQYSHGGQPPEGMNQYYPDGHNDYGY  
 QQPSYPEQGYDRPYEDSSQHYYEGGNSQYGGQQDAYQGPPPPQQGYPPQQQ  
 QYPGQQGYPGQQGYGPSQGGPGPYPNYPQGGQQYGGYRPTQPGPPQP  
 PQRYPYGDQQGYGNYQ

分子量 72 kDa including tags

氨基酸 1 to 418

#### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab116825** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA
	SDS-PAGE
	Western blot

形式 Liquid

补充说明 Protein concentration is above or equal to 0.05 mg/ml.

Best use within three months from the date of receipt of this protein.

## 制备和贮存

### 稳定性和存储

Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.

pH: 8.00

Constituents: 0.3% Glutathione, 0.79% Tris HCl

## 常规信息

### 功能

Appears to function synergistically with RBM14 as a transcriptional coactivator. Isoform 1 and isoform 2 function in nuclear receptor coactivation. Isoform 1 and isoform 2 function in general transcriptional coactivation.

### 组织特异性

Fairly ubiquitously expressed. Expressed in synovial sarcomas and in other human cell lines. The fusion genes SSXT-SSX1 and SSXT-SSX2 are expressed only in synovial sarcomas.

### 疾病相关

Note=A chromosomal aberration involving SS18 may be a cause of synovial sarcoma. Translocation t(X;18)(p11.2;q11.2). The translocation is specifically found in more than 80% of synovial sarcoma. The fusion products SSXT-SSX1 or SSXT-SSX2 are probably responsible for transforming activity. Heterogeneity in the position of the breakpoint can occur (low frequency).


### 序列相似性

Belongs to the SS18 family.

### 细胞定位

Nucleus.

## 图片



12.5% SDS-PAGE stained with Coomassie

Blue showing ab116825 at approximately

72.05 kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish

- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors