

Product datasheet

Recombinant Human STMN2 protein ab89354

1 图像

概述

产品名称	Recombinant人STMN2 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
种属	Human
氨基酸	39 to 179

技术指标

Our [Abpromise guarantee](#) covers the use of **ab89354** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 90 % SDS-PAGE. ab89354 is purified using conventional chromatography techniques.
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 10% Glycerol, 50mM Tris HCl, pH 7.5
--------	--

常规信息

功能	Is a key regulator of neurite extension through regulation of microtubule instability.
组织特异性	Neuron specific.
序列相似性	Belongs to the stathmin family.


翻译后修饰

Sumoylated.

细胞定位

Cytoplasm. Cytoplasm > perinuclear region. Cell projection > growth cone. Membrane. Cell projection > axon. Associated with punctate structures in the perinuclear cytoplasm, axons, and growth cones of developing neurons. SCG10 exists in both soluble and membrane-bound forms.

Recombinant Human STMN2 protein 图像



15% SDS-PAGE showing ab89354 at approximately 16.4kDa (3µg).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors