

Product datasheet

Recombinant Human Tbx3 protein ab114501

1 图像

概述

产品名称	重组人Tbx3蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Wheat germ

氨基酸序列

Accession	O15119
种属	Human
序列	KENGTSDSESSEQAAFNCFAQASSPAASTVGTSNLKDLCPSEGSDAEAE SKEEHGPEACDAAKISTTTSEEPCKRDKGSPAVKAHLFAAERPRDSGRDLK
分子量	37 kDa including tags
氨基酸	331 to 430

技术指标

Our [Abpromise guarantee](#) covers the use of **ab114501** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	ELISA SDS-PAGE Western blot
形式	Liquid
补充说明	Protein concentration is above or equal to 0.05 mg/ml. This protein is best used within three months from the date of receipt.

制备和贮存


稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. pH: 8.00
--------	---

Constituents: 0.3% Glycerol, 0.79% Tris HCl

常规信息

功能	Transcriptional repressor involved in developmental processes. Probably plays a role in limb pattern formation.
组织特异性	Widely expressed.
疾病相关	Defects in TBX3 are the cause of ulnar-mammary syndrome (UMS) [MIM:181450]. UMS is characterized by ulnar ray defects, obesity, hypogenitalism, delayed puberty, hypoplasia of nipples and apocrine glands.
序列相似性	Contains 1 T-box DNA-binding domain.
细胞定位	Nucleus.

图片



ab114501 analysed by 12.5% SDS-PAGE and stained with Coomassie Blue.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors