

Product datasheet

Recombinant Human UBL5 protein ab109343

1 图像

概述

| | |
|------|--------------------------|
| 产品名称 | Recombinant人UBL5 protein |
| 蛋白长度 | Full length protein |

描述

| | |
|----|------------------|
| 性质 | Recombinant |
| 来源 | Escherichia coli |

氨基酸序列

| | |
|-----------|--------------------------------------------------------------------------------------------------------|
| Accession | Q9BZL1 |
| 种属 | Human |
| 序列 | MGSSHHHHHSSGLVPRGSHMIEVVCNDRLGKKVRVKCNTDDTIGDLKKL IAAQTGTRWNKIVLKKWYTIKFDHVSLGDYEIHDGMNLELYQ |
| 分子量 | 11 kDa including tags |
| 氨基酸 | 1 to 73 |
| 标签 | His tag N-Terminus |

技术指标

Our [Abpromise guarantee](#) covers the use of **ab109343** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| | |
|-----|-----------------------------------------------------------------------------|
| 应用 | Mass Spectrometry SDS-PAGE |
| 质谱法 | MALDI-TOF |
| 纯度 | > 90 % SDS-PAGE. ab109343 is purified using conventional chromatography. |
| 形式 | Liquid |

制备和贮存

| | |
|--------|--------------------------------------------------------------------------------------------------|
| 稳定性和存储 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or - |
|--------|--------------------------------------------------------------------------------------------------|

80°C. Avoid freeze / thaw cycle.


Preservative: None

Constituents: 20% Glycerol, 0.1M Sodium chloride, 20mM Tris HCl, 2mM DTT, pH 8.0

常规信息

| | |
|-------|----------------------------------------------------------------------------------------------------------|
| 组织特异性 | Ubiquitous. Highest level of expression in heart, skeletal muscle, kidney, liver, iris and lymphoblasts. |
| 序列相似性 | Contains 1 ubiquitin-like domain. |
| 细胞定位 | Cytoplasm. |

图片



15% SDS-PAGE analysis of 3 µg ab109343.

SDS-PAGE - UBL5 protein (ab109343)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors