

Product datasheet

Recombinant Human beta Defensin 1 protein ab50058

概述

产品名称	Recombinant人beta Defensin 1 protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

种属	Human
----	-------

技术指标

Our [Abpromise guarantee](#) covers the use of **ab50058** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
----	--------------------------------

内毒素水平	< 0.100 Eu/μg
-------	---------------

纯度	> 95 % SDS-PAGE.
----	------------------

形式	Lyophilised
----	-------------

制备和贮存

稳定性和存储 Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.

Preservative: None

Endotoxin level is less than 0.1 ng per μg (1EU/μg).

复溶	Reconstitute to 1mg/ml using sterile water.
----	---

常规信息

功能	Has bactericidal activity.
----	----------------------------

组织特异性	Plasma.
序列相似性	Belongs to the beta-defensin family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors