

Product datasheet

Recombinant KIR2DL1 protein ab51212

1 图像

概述

产品名称	Recombinant KIR2DL1 protein
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli

技术指标

Our [Abpromise guarantee](#) covers the use of **ab51212** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 95 % SDS-PAGE. ab51212 was purified by FPLC gel-filtration chromatography, after refolding of the isolated inclusion bodies in a redox buffer.
形式	Liquid
补充说明	ab51212 corresponds to the extracellular domain of KIR2DL1 and was cloned into an E. coli expression vector and overexpressed as insoluble protein aggregates (inclusion bodies).

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 20mM Tris HCl, pH 7.5
--------	--

常规信息

功能	Receptor on natural killer (NK) cells for HLA-C alleles. Inhibits the activity of NK cells thus preventing cell lysis.
----	--


序列相似性

Belongs to the immunoglobulin superfamily.
Contains 2 Ig-like C2-type (immunoglobulin-like) domains.

细胞定位

Cell membrane.

图片



14% SDS-PAGE gel loaded with ab51212,
recombinant human KIR2DL1 protein
(extracellular domain).

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors