

Recombinant Mouse DR6 protein ab109063

描述

产品名称	重组小鼠DR6蛋白
纯度	> 95 % SDS-PAGE.
内毒素水平	< 0.100 Eu/μg
表达系统	HEK 293 cells
Accession	<u>Q9EPU5</u>
蛋白长度	Protein fragment
无动物成分	No
性质	Recombinant
种属	Mouse
预测分子量	80 kDa
氨基酸	42 to 348

技术指标

Our **Abpromise guarantee** covers the use of **ab109063** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
形式	Lyophilized
补充说明	PBS containing at least 0.1% BSA should be used for further dilutions.

制备和贮存

稳定性和存储	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C. Constituent: PBS
复溶	Reconstitute with 50 μl sterile water.

常规信息

功能	May activate NF-kappa-B and promote apoptosis. May activate JNK and be involved in T-cell
----	---

differentiation. Required for both normal cell body death and axonal pruning. Trophic-factor deprivation triggers the cleavage of surface APP by beta-secretase to release sAPP-beta which is further cleaved to release an N-terminal fragment of APP (N-APP). N-APP binds TNFRSF21 triggering caspase activation and degeneration of both neuronal cell bodies (via caspase-3) and axons (via caspase-6).

组织特异性

Highly expressed in heart, brain, placenta, pancreas, lymph node, thymus and prostate. Detected at lower levels in lung, skeletal muscle, kidney, testis, uterus, small intestine, colon, spleen, bone marrow and fetal liver. Very low levels were found in adult liver and peripheral blood leukocytes.

序列相似性

Contains 1 death domain.
Contains 4 TNFR-Cys repeats.

细胞定位

Membrane.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors