

Product datasheet

Recombinant Mouse IL-17E protein ab119730

概述

产品名称	重组小鼠IL-17E蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	Q8VHH8
种属	Mouse
序列	VSLRIQEGCSHLPSCCPSKEQPPEEWLKWSSASVSPPEP LSHTTHAESC RASKDGPLNSRAISPWSYELDRDLNRVPQDLYHARCLCPH CVSLQTGSHM DPLGNSVPLYHNQTVFYRRPCHGEEGTRRYCLERRLYRVS LACVCVRPRV MA
分子量	35 kDa
氨基酸	17 to 169

技术指标

Our [Abpromise guarantee](#) covers the use of **ab119730** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	HPLC Functional Studies
纯度	> 95 % SDS-PAGE. Purity by SDS-PAGE and HPLC: = or > 95% Endotoxin level is < 0.01 ng per µg of ab119730.
形式	Lyophilised
补充说明	The lyophilized ab119730 is best-stored desiccated at -20 °C. Upon reconstitution the product should be aliquoted and stored at -20 °C. It is recommended to add a carrier protein (0.1% HSA or BSA) for long term storage.

制备和贮存

稳定性和存储	Shipped at 4°C. Store at -20°C.
复溶	Centrifuge vial before opening. When reconstituting the product, gently pipet and wash down the sides of the vial to ensure full recovery of the protein into solution. It is recommended to reconstitute the lyophilized product with sterile 5 mM Acetic acid to a concentration of 0.1 mg/ml, which can be further diluted into other aqueous solutions.
常规信息	
功能	Induces activation of NF-kappa-B and stimulates production of the proinflammatory chemokine IL-8. Proinflammatory cytokine favoring Th2-type immune responses.
组织特异性	Expressed at low levels in several tissues, including brain, kidney, lung, prostate, testis, spinal cord, adrenal gland, and trachea.
序列相似性	Belongs to the IL-17 family.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors