

## Product datasheet

# Recombinant Mouse IL9 protein ab51918

### 概述

---

产品名称	重组小鼠IL9蛋白
蛋白长度	Full length protein

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

种属	Mouse
----	-------

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab51918** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
纯度	> 95 % SDS-PAGE. >98% by SDS-PAGE and HPLC analyses.
形式	Lyophilised
补充说明	The reconstituted solution can be diluted into other aqueous buffers and stored at 4°C for 1 week or -20°C for future use.

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Endotoxin level is <0.1 ng per µg.
复溶	Reconstitute in sterile distilled water to a concentration of 0.1-1.0 mg/ml.

### 常规信息

---

功能	Supports IL-2 independent and IL-4 independent growth of helper T-cells.
序列相似性	Belongs to the IL-7/IL-9 family.
细胞定位	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors