

Product datasheet

Recombinant Mouse Resistin protein ab50178

概述

产品名称	Recombinant小鼠Resistin protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
种属	Mouse
序列	SSMPLCPIDE AIDKKIKQDF NSLFPNAIKN IGLNCWTVSS RGKLASCPEG TAVLSCSCGS ACGSWDIREE KVCHCQCARI DWTAARCCKL QVAS

技术指标

Our [Abpromise guarantee](#) covers the use of **ab50178** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
内毒素水平	< 0.100 Eu/μg
形式	Lyophilised

制备和贮存

稳定性和存储	Shipped at 4°C. The lyophilized protein is stable for a few weeks at room temperature. Store at -20°C long term.
复溶	Reconstitute to 1mg/ml using sterile water.

常规信息

功能	Hormone that seems to suppress insulin ability to stimulate glucose uptake into adipose cells. Potentially links obesity to diabetes.
组织特异性	Expressed only in fatty tissues.
序列相似性	Belongs to the resistin/FIZZ family.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors