

Product datasheet

Recombinant Mouse TIE1 protein ab54433

概述

产品名称	Recombinant小鼠TIE1 protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	CHO cells

氨基酸序列

种属	Mouse
----	-------

技术指标

Our [Abpromise guarantee](#) covers the use of **ab54433** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE
纯度	> 90 % SDS-PAGE. Purity: > 90%, by SDS-PAGE and visualised by silver stain. Endotoxin level: < 0.1 ng per µg sTIE-1/Fc
形式	Lyophilised

制备和贮存

稳定性和存储	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. Preservative: None Constituents: PBS
复溶	The protein can be reconstituted in ddH ₂ O or PBS to a concentration of 100µg/ml.

常规信息

功能	Transmembrane tyrosine-protein kinase that may modulate TEK/TIE2 activity and contribute to the regulation of angiogenesis.
----	---

组织特异性	Specifically expressed in developing vascular endothelial cells.
序列相似性	Belongs to the protein kinase superfamily. Tyr protein kinase family. Tie subfamily. Contains 3 EGF-like domains. Contains 3 fibronectin type-III domains. Contains 2 Ig-like C2-type (immunoglobulin-like) domains. Contains 1 protein kinase domain.
翻译后修饰	Phosphorylated on tyrosine residues in response to ANGPT1, most likely by TEK/TIE2.
细胞定位	Cell membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors