

## Product datasheet

# Recombinant Rat Choline Acetyltransferase protein ab119467

### 概述

---

产品名称	重组大鼠Choline Acetyltransferase蛋白
蛋白长度	Full length protein

### 描述

---

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

Accession	<a href="#">P32738</a>
种属	Rat
分子量	72 kDa
氨基酸	1 to 640

### 技术指标

---

Our [Abpromise guarantee](#) covers the use of **ab119467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Immunohistochemistry (Frozen sections) Blocking
----	--

形式	Liquid
----	--------

### 制备和贮存

---

稳定性和存储	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C. Constituents: 89% PBS, 10% Glycerol
--------	---

### 常规信息

---

功能	Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.
----	--

疾病相关	Myasthenic syndrome, congenital, 6, presynaptic
------	---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors