

Product datasheet

Recombinant Rat Choline Acetyltransferase protein ab119467

概述

产品名称	Recombinant大鼠Choline Acetyltransferase protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

Accession	P32738
种属	Rat
分子量	72 kDa
氨基酸	1 to 640

技术指标

Our [Abpromise guarantee](#) covers the use of **ab119467** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Immunohistochemistry (Frozen sections) Blocking
----	----------------------------------------------------

形式	Liquid
----	--------

制备和贮存

稳定性和存储	Shipped at 4°C. Store at 4°C (up to 6 months). Store at -20°C. Constituents: 89% PBS, 10% Glycerol
--------	-------------------------------------------------------------------------------------------------------

常规信息

功能	Catalyzes the reversible synthesis of acetylcholine (ACh) from acetyl CoA and choline at cholinergic synapses.
----	----------------------------------------------------------------------------------------------------------------

疾病相关	Myasthenic syndrome, congenital, 6, presynaptic
------	-------------------------------------------------

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors