

Product datasheet

Recombinant human Endostatin protein ab73653

概述

产品名称	Recombinant人Endostatin protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Pichia pastoris
氨基酸序列	
种属	Human
序列	Amino Acid Composition is in total agreement with the expected amino acid composition of native human Endostatin.

技术指标

Our [Abpromise guarantee](#) covers the use of **ab73653** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生物活性	The activity calculated by ECE migration inhibition was found to be 50,000IU/mg.
应用	SDS-PAGE
纯度	> 95 % SDS-PAGE. ab73653 is purified by proprietary chromatographic techniques. Purity is greater than 98.0% as determined by RP-HPLC and SDS-PAGE.
形式	Liquid

制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

常规信息

功能	COLA18A probably plays a major role in determining the retinal structure as well as in the closure of the neural tube. Endostatin potently inhibits endothelial cell proliferation and angiogenesis. May inhibit angiogenesis by binding to the heparan sulfate proteoglycans involved in growth factor signaling.
组织特异性	Present in multiple organs with highest levels in liver, lung and kidney.
疾病相关	Defects in COL18A1 are a cause of Knobloch syndrome (KNO) [MIM:267750]. KNO is an autosomal recessive disorder defined by the occurrence of high myopia, vitreoretinal degeneration with retinal detachment, macular abnormalities and occipital encephalocele.
序列相似性	Belongs to the multiplexin collagen family. Contains 1 FZ (frizzled) domain. Contains 1 TSP N-terminal (TSPN) domain.
翻译后修饰	Prolines at the third position of the tripeptide repeating unit (G-X-Y) are hydroxylated in some or all of the chains.
细胞定位	Secreted > extracellular space > extracellular matrix.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors