

Product datasheet

Recombinant human PTPRS protein ab61149

2 图像

概述

产品名称	重组人PTPRS蛋白
蛋白长度	Protein fragment

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
种属	Human
分子量	63 kDa including tags
氨基酸	883 to 1210
标签	GST tag N-Terminus
额外的序列信息	GeneBank Accession number: NM_130855

技术指标

Our [Abpromise guarantee](#) covers the use of **ab61149** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	SDS-PAGE Functional Studies
形式	Liquid


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 30% Glycerol, 50mM Sodium chloride, 20mM MOPS, 0.25mM DTT, 0.1mM PMSF, pH 7.5 This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	--

## 常规信息


功能	Interacts with LAR-interacting protein LIP.1.
组织特异性	Detected in all tissues tested except for placenta and liver.
序列相似性	Belongs to the protein-tyrosine phosphatase family. Receptor class 2A subfamily. Contains 8 fibronectin type-III domains. Contains 3 Ig-like C2-type (immunoglobulin-like) domains. Contains 2 tyrosine-protein phosphatase domains.
翻译后修饰	A cleavage occurs, separating the extracellular domain from the transmembrane segment. This process called 'ectodomain shedding' is thought to be involved in receptor desensitization, signal transduction and/or membrane localization.
细胞定位	Membrane.

## 图片



Sample Kinase Activity Plot.

Functional Studies - PTPRS protein (Tagged)  
(ab61149)



ab61149 on SDS-PAGE, MW ~63 kDa.

SDS-PAGE - PTPRS protein (Tagged) (ab61149)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors