

Product datasheet

Recombinant human Prolactin protein ab119163

概述

产品名称	Recombinant人Prolactin protein
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli
氨基酸序列	
Accession	P01236
种属	Human
分子量	23 kDa
氨基酸	29 to 227

技术指标

Our [Abpromise guarantee](#) covers the use of **ab119163** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生物活性	The ED50 is < 0.065ng/ml, as determined by the dose dependant stimulation of the proliferation of Rat lymphoma Nb211.
应用	Functional Studies SDS-PAGE HPLC
纯度	> 95 % SDS-PAGE.
形式	Lyophilised

制备和贮存

稳定性和存储	Shipped at 4°C. After reconstitution store at -20°C. Avoid freeze / thaw cycles. Constituent: 99% PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute with 100 ul distilled water.

常规信息

相关性	Various hormones are secreted from the anterior pituitary during development and growth. Prolactin is a growth factor secreted by the anterior pituitary that is necessary for the proliferation and differentiation of the mammary glands.
细胞定位	Secreted

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors