

Product datasheet

Recombinant human TSSK2 protein ab70075

2 图像

概述

产品名称	重组人TSSK2蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Baculovirus infected Sf9 cells
氨基酸序列	
种属	Human

技术指标

Our [Abpromise guarantee](#) covers the use of **ab70075** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
纯度	> 90 % SDS-PAGE.
形式	Liquid
补充说明	<a href="#">ab64311</a> (Myelin Basic Protein protein) can be utilized as a substrate for assessing kinase activity


制备和贮存

稳定性和存储	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. Preservative: None Constituents: 25% Glycerol, 50mM Tris HCl, 150mM Sodium chloride, 0.25mM DTT, 0.1mM EGTA, 0.1mM EDTA, 0.1mM PMSF, pH 7.5 This product is an active protein and may elicit a biological response in vivo, handle with caution.
--------	---


## 常规信息

<b>功能</b>	Testis-specific serine/threonine-protein kinase required during spermatid development. Phosphorylates TSKS at 'Ser-288' and SPAG16. Involved in the late stages of spermatogenesis, during the reconstruction of the cytoplasm. During spermatogenesis, required for the transformation of a ring-shaped structure around the base of the flagellum originating from the chromatoid body.
<b>组织特异性</b>	Testis-specific. Present in mature spermatozoa (at protein level).
<b>序列相似性</b>	Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. Contains 1 protein kinase domain.
<b>翻译后修饰</b>	Autophosphorylated.
<b>细胞定位</b>	Cytoplasm. Cytoplasm > cytoskeleton > centrosome > centriole. Present in the cytoplasm of elongating spermatids. In spermatozoa, localizes in the equatorial segment, neck, the midpiece and in a specific sperm head compartment (By similarity). In spermatids, concentrates in centrioles during flagellogenesis.

## 图片



Kinase activity assay of ab70075. The specific activity of TSSK2 was determined to be 136 nmol/min/mg.



SDS-PAGE of ab70075. Molecular weight 66kDa.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors