

## Product datasheet

# Recombinant mouse IL19 protein ab51987

### 概述

产品名称	Recombinant小鼠IL19 protein
蛋白长度	Full length protein

### 描述

性质	Recombinant
来源	Escherichia coli

### 氨基酸序列

种属	Mouse
----	-------

### 技术指标

Our [Abpromise guarantee](#) covers the use of **ab51987** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Functional Studies SDS-PAGE
纯度	> 95 % SDS-PAGE. >95% by SDS-PAGE and HPLC analyses
形式	Lyophilised
补充说明	The reconstituted solution can be diluted into aqueous buffers and stored at 4°C for 1 week or –20°C for future use.

### 制备和贮存

稳定性和存储	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. Preservative: None Endotoxin level is <0.1 ng per µg of IL19. This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	Reconstitute in water to a concentration of 0.1-0.5 mg/ml.

### 常规信息

---

功能	May play some important roles in inflammatory responses. Up-regulates IL-6 and TNF-alpha and induces apoptosis.
序列相似性	Belongs to the IL-10 family.
细胞定位	Secreted.

---

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

---

### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

---

### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors