

Product datasheet

Recombinant rat GRO gamma protein ab125819

概述

产品名称	重组大鼠GRO gamma蛋白
蛋白长度	Full length protein

描述

性质	Recombinant
来源	Escherichia coli

氨基酸序列

Accession [Q10746](#)

种属 Rat

序列 RELRCQCLKT LPRVDFENIQ SLTVTPPGPH
CTQTEVIATL KDGQEVCLNP QAPRLQKIIQ
KLLKSDKSS

分子量 8 kDa

氨基酸 33 to 101

额外的序列信息 GRO-g Rat Recombinant produced in E.Coli is a single,non-glycosylated, polypeptide chain containing 69 amino acids and having a molecular mass of 7.8kDa.

技术指标

Our [Abpromise guarantee](#) covers the use of **ab125819** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

生物活性 Determined by its ability to chemoattract 293 transfected CXCR2 cells using a concentration range of 10-100 ng/ml, corresponding to a specific activity of 10,000-100,000 units/mg.

应用 Functional Studies
SDS-PAGE
HPLC

纯度 > 95 % SDS-PAGE.
ab125819 is purified by proprietary chromatographic techniques. Purity is greater than 97.0% as determined by RP-HPLC and SDS-PAGE.

形式 Lyophilised

补充说明	Upon reconstitution ab125819 should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.
制备和贮存	
稳定性和存储	Shipped at 4°C. Store at -20°C. Store under desiccating conditions. pH: 7.40 Constituent: 99% PBS This product is an active protein and may elicit a biological response in vivo, handle with caution.
复溶	It is recommended to reconstitute the lyophilized ab125819 in sterile 18Mcm H ₂ O not less than 100µg/ml, which can then be further diluted to other aqueous solutions.
常规信息	
功能	Ligand for CXCR2 (By similarity). Has chemotactic activity for neutrophils. May play a role in inflammation and exert its effects on endothelial cells in an autocrine fashion. In vitro, the processed form GRO-gamma(5-73) shows a fivefold higher chemotactic activity for neutrophilic granulocytes.
序列相似性	Belongs to the intercrine alpha (chemokine CxC) family.
翻译后修饰	N-terminal processed form GRO-gamma(5-73) is produced by proteolytic cleavage after secretion from peripheral blood monocytes.
细胞定位	Secreted.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors