


Product datasheet

Anti-Rhophilin 2 antibody ab79011

1 图像

概述

产品名称	Anti-Rhophilin 2抗体
描述	山羊多克隆抗体to Rhophilin 2
经测试应用	适用于: ELISA, IHC-P
种属反应性	与反应: Mouse, Rat, Human 预测可用于: Cow 
免疫原	Synthetic peptide: C-TNKNRQKSAST , corresponding to internal sequence amino acids 645-655 of Human Rhophilin 2 (NP_149094.2). Run BLAST with Run BLAST with
阳性对照	Kidney tissue.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 0.5% BSA, Tris saline, pH 7.3
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab79011** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ELISA		1/32000.
IHC-P		Use a concentration of 5 µg/ml.

靶标

功能	Binds specifically to GTP-Rho. May function in a Rho pathway to limit stress fiber formation and/or increase the turnover of F-actin structures in the absence of high levels of RhoA activity.
组织特异性	Widely expressed. Highly expressed in prostate, trachea, stomach, colon, thyroid and pancreas. Expressed at lower level in brain, spinal cord, kidney, placenta and liver.
序列相似性	Belongs to the RHPN family. Contains 1 BRO1 domain. Contains 1 PDZ (DHR) domain. Contains 1 REM (Hr1) repeat.
细胞定位	Cytoplasm > perinuclear region.

Anti-Rhopilin 2 antibody 图像



ab79011, at 5µg/ml, staining Rhophilin 2 in kidney tissue by Immunohistochemistry using formalin-fixed, paraffin-embedded tissue.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Rhophilin 2 antibody (ab79011)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors