

Product datasheet

Anti-SCAND1 antibody ab64828

2 图像

概述

产品名称	Anti-SCAND1抗体
描述	兔多克隆抗体to SCAND1
宿主	Rabbit
特异性	ab64828 detects endogenous levels of total SCAND1 protein.
经测试应用	适用于: ELISA, WB, ICC/IF
种属反应性	与反应: Human
免疫原	Synthetic peptide derived from an internal sequence of human SCAND1.
阳性对照	Jurkat cell extracts and HepG2 cells.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at -20°C. Stable for 12 months at -20°C.
存储溶液	Preservative: 0.02% Sodium Azide Constituents: 50% Glycerol, PBS (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), 150mM Sodium chloride, pH 7.4
纯度	Immunogen affinity purified
纯化说明	ab64828 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab64828** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
ELISA		1/5000.

应用	Ab评论	说明
WB		1/500 - 1/1000. Predicted molecular weight: 19 kDa.
ICC/IF		1/500 - 1/1000.


## 靶标

功能	May regulate transcriptional activity.
序列相似性	Contains 1 SCAN box domain.
细胞定位	Nucleus.

## 图片



ab64828, at 1/500 dilution, staining SCAND1 in HepG2 cells by Immunofluorescence, in the absence (left image) or presence (right image) of the immunising peptide.



**All lanes :** Anti-SCAND1 antibody (ab64828) at 1/500 dilution

**Lane 1 :** Jurkat cell extract at 5 µg

**Lane 2 :** Jurkat cell extract plus immunising peptide at 5 µg

**Predicted band size:** 19 kDa

**Observed band size:** 25 kDa

**Additional bands at:** 70 kDa (possible isoform)

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors