


Product datasheet

Anti-SCRN2 antibody ab117731

2 图像

概述

产品名称	Anti-SCRN2抗体
描述	兔多克隆抗体to SCRN2
宿主	Rabbit
经测试应用	适用于: WB, ELISA, IHC-P, ICC
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide conjugated to KLH, corresponding to a region from the C-terminus of Human SCRN2.
阳性对照	Human Spleen and Human Thyroid tissues

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.02% Sodium azide Constituent: 99% PBS
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab117731** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.5 - 1 µg/ml. Predicted molecular weight: 46 kDa.
ELISA		Use at an assay dependent concentration.


应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.
ICC		Use at an assay dependent concentration.

靶标

序列相似性


Belongs to the peptidase C69 family. Secernin subfamily.

图片



ab117731 at 5 µg/ml staining SCR3 in Formalin-Fixed, Paraffin-Embedded Human Spleen.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCR2 antibody (ab117731)



ab117731 at 5 µg/ml staining SCR3 in Formalin-Fixed, Paraffin-Embedded Human Thyroid.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SCR2 antibody (ab117731)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors