


Product datasheet

Anti-SGK3 antibody ab61818

1 图像

概述

产品名称	Anti-SGK3抗体
描述	兔多克隆抗体to SGK3
宿主	Rabbit
经测试应用	适用于: IHC-P
种属反应性	与反应: Human 预测可用于: Mouse, Rat 
免疫原	Synthetic peptide corresponding to Human SGK3 (N terminal) conjugated to Keyhole Limpet Haemocyanin (KLH).
阳性对照	Human Pancreas, Islet tissue
常规说明	For short term storage: at -4°C

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Avoid freeze / thaw cycle.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, pH 7.7
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab61818** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
IHC-P		Use a concentration of 5 µg/ml.

靶标

功能	Involved in the activation of potassium channels. Mediates cell IL-3-dependent survival signals.
组织特异性	Expressed in most tissues with highest levels in pancreas, kidney liver, heart and brain and lower levels in lung, placenta and skeletal muscle.
序列相似性	Belongs to the protein kinase superfamily. AGC Ser/Thr protein kinase family. Contains 1 AGC-kinase C-terminal domain. Contains 1 protein kinase domain. Contains 1 PX (phox homology) domain.
细胞定位	Cytoplasmic vesicle. Early endosome. Localized in vesicle-like structures and in the early endosome.

图片



Immunohistochemistry analysis of human Pancreas (Formalin-Fixed Paraffin-Embedded) using ab61818 at 5µg/ml.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-SGK3 antibody (ab61818)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors