

Product datasheet

Anti-SLC41A3 antibody ab80462

1 图像

概述

产品名称	Anti-SLC41A3抗体
描述	兔多克隆抗体to SLC41A3
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Caenorhabditis elegans, Drosophila melanogaster, Zebrafish
免疫原	Synthetic peptide corresponding to Human SLC41A3 aa 434-483 (C terminal). Sequence: WHQALDPDNHCIPYLTGLGDLLGSSSVGHATAAVPRRCTASPGWGLIQPFI Run BLAST with Run BLAST with
阳性对照	HeLa cell lysate.

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Preservative: None Constituents: 2% Sucrose, PBS
纯度	Immunogen affinity purified
纯化说明	ab80462 is purified by a peptide affinity chromatography method.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab80462** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab 评论	说明
----	-------	----

WB Use a concentration of 1 µg/ml. Detects a band of approximately 55 kDa (predicted molecular weight: 55 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


ELISA Use at an assay dependent concentration.
ELISA titre using peptide based assay: 1/62500.

靶标

相关性 SLC41A3 belongs to the SLC41A transporter family. The exact function of SLC41A3 remains unknown.

细胞定位 Cell membrane; Multi-pass membrane protein.

图片



Anti-SLC41A3 antibody (ab80462) at 1 µg/ml
+ HeLa cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

Predicted band size: 55 kDa

Observed band size: 55 kDa

Additional bands at: 35 kDa, 40 kDa, 85 kDa. We are unsure as to the identity of these extra bands.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors