

Product datasheet

Anti-SLIRP antibody ab51523

★★★★★ 1 Abreviews 6 References 1 图像

概述

| | |
|-------|---|
| 产品名称 | Anti-SLIRP抗体 |
| 描述 | 兔多克隆抗体to SLIRP |
| 宿主 | Rabbit |
| 经测试应用 | 适用于: IHC-Fr, WB |
| 种属反应性 | 与反应: Human 不与反应: Mouse |
| 免疫原 | Synthetic peptide made to a C terminal region of human SLIRP. |
| 阳性对照 | Human liver lysate |

性能

| | |
|------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. |
| 存储溶液 | Preservative: 0.05% Sodium Azide Constituents: Tris glycine, 150mM Sodium chloride |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab51523** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|--------|-------|--|
| IHC-Fr | ★★★★★ | Use at an assay dependent concentration. |


| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use a concentration of 0.5 - 2 µg/ml. Detects a band of approximately 12 kDa (predicted molecular weight: 12 kDa). |

靶标

相关性 SLIRP is a RNA binding protein that acts as a nuclear receptor corepressor. It may act by binding the SRA RNA, and repressing the SRA-mediated nuclear receptor coactivation. SLIRP binds the STR7 loop of SRA RNA. It is also able to repress glucocorticoid (GR), androgen (AR), thyroid (TR) and VDR-mediated transactivation. It contains 1 RRM (RNA recognition motif) domain. SLIRP is predominantly found in mitochondria with a limited role in the nucleus where it is recruited to nuclear receptor target promoters.

细胞定位 Mitochondrial and Nuclear; Predominantly mitochondrial.

图片



Anti-SLIRP antibody (ab51523) at 0.6 µg/ml + human liver lysate

Predicted band size: 12 kDa

Observed band size: 12 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors