

## Product datasheet

# Anti-SNRPC antibody [4H12] ab122901

★★★★☆ 1 Abreviews 1 References 3 图像

### 概述

产品名称	Anti-SNRPC抗体[4H12]
描述	大鼠单克隆抗体[4H12] to SNRPC
宿主	Rat
经测试应用	适用于: WB, IP, ELISA, ICC/IF
种属反应性	与反应: Mouse, Human, Monkey
免疫原	Recombinant full length Mouse SNRPC with a proprietary tag.
阳性对照	L929 cells; HeLa total cell extract.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: 49% PBS, 50% Glycerol
纯化说明	ab122901 was produced from the hybridoma cultured in serum-free medium and purified under mild conditions by a propriety chromatography processes and filter sterilized.
克隆	单克隆
克隆编号	4H12
同种型	IgG2a
轻链类型	kappa

### 应用

Our [Abpromise guarantee](#) covers the use of **ab122901** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


应用	Ab评论	说明
WB	★★★★☆	Use at an assay dependent dilution. Predicted molecular weight: 17 kDa.

应用	Ab评论	说明
IP		Use at an assay dependent dilution.
ELISA		Use at an assay dependent dilution.
ICC/IF		Use at an assay dependent dilution.


## 靶标

功能	This protein is associated with snRNP U1.
序列相似性	Belongs to the U1 small nuclear ribonucleoprotein C family. Contains 1 matrin-type zinc finger.
细胞定位	Nucleus.

## 图片




ab122901 staining SNRPC in Mouse L929 cells by Immunofluorescence (green). Right image shows the cells DNA stained with Hoechst.



Anti-SNRPC antibody [4H12] (ab122901) + HeLa total cell extract.

**Predicted band size: 17 kDa**



Immunoprecipitation - Anti-SNRPC antibody [4H12]  
(ab122901)

Immunoprecipitation of SNRPC with  
ab122901.

Lane 1: no extract.

Lane 2: HeLa total cell extract.

HL and LC: heavy and light chains.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors