

## Product datasheet

# Anti-SPNS antibody ab59978

### 概述

产品名称	Anti-SPNS抗体
描述	兔多克隆抗体to SPNS
宿主	Rabbit
特异性	This antibody reacts specifically with SPNS proteins. It cross reacts with the mouse SPNS1, SPNS2 and SPNS3 and human SPNS1 isoform 1, 2, 3 and 4.
经测试应用	适用于: WB, IHC-P, IHC-Fr
种属反应性	与反应: Mouse, Human
免疫原	Synthetic peptide derived from mouse common domain of SPNS proteins.

### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term.
存储溶液	Preservative: None Constituents: Whole serum
纯度	Whole antiserum
克隆	多克隆
同种型	IgG

### 应用

Our [Abpromise guarantee](#) covers the use of **ab59978** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/2000. Predicted molecular weight: 58 kDa.
IHC-P		1/100 - 1/500.
IHC-Fr		1/100 - 1/500.

## 靶标

---

**相关性** Sphingolipid transporter required for migration of myocardial precursors. Transports sphingosine 1-phosphate (S1P), a secreted lipid mediator that plays critical roles in cardiovascular, immunological, and neural development and function. Mediates the export of S1P from cells in the extraembryonic yolk syncytial layer (YSL), thereby regulating myocardial precursor migration.

**细胞定位** Membrane; Multi pass membrane protein

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors