

Product datasheet

Anti-STARD4 antibody ab111402

1 图像

概述

| | |
|--------------|---|
| 产品名称 | Anti-STARD4抗体 |
| 描述 | 兔多克隆抗体to STARD4 |
| 经测试应用 | 适用于: WB |
| 种属反应性 | 与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish |
| 免疫原 | Synthetic peptide corresponding to a region within internal sequence amino acids 96-145 (CCVMRYTTAG QLWNIISPRE FVDFSYTVGY KEGLLSGIS LDWDEKRPEF) of Human STARD4 (NP_631903). Run BLAST with ExPASy Run BLAST with NCBI |
| 阳性对照 | HepG2 cell lysate. |

性能

| | |
|-------------|--|
| 形式 | Liquid |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| 存储溶液 | Constituents: PBS, 2% Sucrose |
| 纯度 | Immunogen affinity purified |
| 克隆 | 多克隆 |
| 同种型 | IgG |

应用

Our [Abpromise guarantee](#) covers the use of **ab111402** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明 |
|----|------|--|
| WB | | Use a concentration of 1 µg/ml. Predicted molecular weight: 24 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |

靶标

| | |
|-------|---|
| 功能 | May be involved in the intracellular transport of sterols or other lipids. May bind cholesterol or other sterols. |
| 序列相似性 | Contains 1 START domain. |

图片



Anti-STARD4 antibody (ab111402) at 1 μ g/ml
+ HepG2 cell lysate at 10 μ g

Predicted band size : 24 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors