

Product datasheet

Anti-STARD4 antibody ab111402

1 图像

概述

产品名称	Anti-STARD4抗体
描述	兔多克隆抗体to STARD4
宿主	Rabbit
经测试应用	适用于: WB
种属反应性	与反应: Human 预测可用于: Mouse, Rat, Rabbit, Horse, Chicken, Guinea pig, Cow, Cat, Dog, Pig, Zebrafish



免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 96-145 (CCVMRYTTAG QLWNIISPRE FVDFSYTVGY KEGLLSGGIS LDWDEKRPEF) of Human STARD4 (NP_631903).
-----	-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

Run BLAST with [ExPASy](#)
 Run BLAST with [NCBI](#)

阳性对照	HepG2 cell lysate.
------	--------------------

性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
存储溶液	Constituents: PBS, 2% Sucrose
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab111402** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明


应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 24 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

靶标

功能 May be involved in the intracellular transport of sterols or other lipids. May bind cholesterol or other sterols.

序列相似性 Contains 1 START domain.

图片



Anti-STARD4 antibody (ab111402) at 1 µg/ml
+ HepG2 cell lysate at 10 µg

Predicted band size: 24 kDa

Gel concentration: 12%.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors