

Product datasheet

Anti-STRBP antibody ab111567

1 图像

概述

|       |  |
|-------|--|
| 产品名称  | Anti-STRBP抗体   |
| 描述    | 兔多克隆抗体to STRBP   |
| 宿主    | Rabbit   |
| 经测试应用 | 适用于: WB  |
| 种属反应性 | 与反应: Human<br>预测可用于: Mouse, Rat, Cow  |
| 免疫原   | Recombinant fragment, corresponding to a region within amino acids 1-165 of Human STRBP.                               |
| 阳性对照  | NT2D1, IMR32 and MCF7 whole cell lysates.  |

性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.           |
| 存储溶液 | pH: 7.00<br>Preservative: 0.01% Thimerosal (merthiolate)<br>Constituents: 20% Glycerol, 1.21% Tris, 0.75% Glycine |
| 纯度   | Immunogen affinity purified   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

应用

Our [Abpromise guarantee](#) covers the use of **ab111567** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用 | Ab评论 | 说明  |
|----|------|---|
| WB |      | 1/500 - 1/3000. Predicted molecular weight: 74 kDa. |

## 靶标

|       |  |
|-------|--|
| 功能    | Involves in spermatogenesis and sperm function. Plays a role in regulation of cell growth. Binds to double-stranded DNA and RNA. Binds most efficiently to poly(I:C) RNA than to poly(dI:dC) DNA. Binds also to single-stranded poly(G) RNA. Binds non-specifically to the mRNA PRM1 3'-UTR and adenovirus VA RNA. |
| 序列相似性 | Contains 2 DRBM (double-stranded RNA-binding) domains.<br>Contains 1 DZF domain.   |
| 细胞定位  | Cytoplasm. Microtubule-associated that localizes to the manchette in developing spermatids.  |

## 图片



Anti-STRBP antibody (ab111567) at 1/1000 dilution + IMR32 whole cell lysate at 30 µg/ml

**Predicted band size:** 74 kDa

7.5% SDS-PAGE.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors