

Product datasheet

Schneider L2 whole cell lysate ab14893

概述

| | |
|------|--|
| 产品名称 | Schneider L2全细胞裂解液 |
| 常规说明 | Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA). To ensure stability, we recommend aliquoting the extracts into single-use fractions and then storing them at -80°C. |

性能

| | |
|------|---|
| 形式 | Liquid |
| 存放说明 | Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles. |
| 存储溶液 | pH: 7.90 Constituents: 0.48% HEPES, 0.75% Potassium chloride, 0.01% Magnesium chloride, 20% Glycerol, 0.009% PMSF, 0.008% DTT |
| 纯度 | Whole Cell Lysate |
| 背景 | This permanent line is derived from a primary culture of late stage (L2) (20-24 hours old) <i>Drosophila melanogaster</i> embryos (Schneider, 1972) prepared in Scheiders insect medium. It is a positive control for the transcription factor GAGA ab11962 . |

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors