

Product datasheet

Schneider L2 whole cell lysate ab14893

概述

产品名称	Schneider L2全细胞裂解液
常规说明	Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA). To ensure stability, we recommend aliquoting the extracts into single-use fractions and then storing them at -80°C.

性能

形式	Liquid
存放说明	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
存储溶液	pH: 7.90 Constituents: 0.48% HEPES, 0.75% Potassium chloride, 0.01% Magnesium chloride, 20% Glycerol, 0.009% PMSF, 0.008% DTT
纯度	Whole Cell Lysate
背景	This permanent line is derived from a primary culture of late stage (L2) (20-24 hours old) <i>Drosophila melanogaster</i> embryos (Schneider, 1972) prepared in Scheiders insect medium. It is a positive control for the transcription factor GAGA ab11962 .

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors