

Product datasheet

Anti-SerpinB2 antibody ab105739

1 图像

概述

产品名称	Anti-SerpinB2抗体
描述	兔多克隆抗体to SerpinB2
宿主	Rabbit
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Human 预测可用于: Monkey 
免疫原	Synthetic peptide corresponding to internal sequence amino acids of Human SerpinB2 (NP_002566.1).
阳性对照	HeLa whole cell lysates.
常规说明	Expiration date is one (1) year from date of opening product.

性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
存储溶液	pH: 7.20 Preservative: 0.01% Sodium azide Constituents: 0.42% Potassium phosphate, 0.88% Sodium chloride
纯度	Protein A purified
纯化说明	ab105739 is affinity purified from monospecific antiserum by immunoaffinity purification.
克隆	多克隆
同种型	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab105739** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		1/500 - 1/2000. Detects a band of approximately 42 kDa (predicted molecular weight: 47 kDa).
ELISA		1/100000 - 1/150000.

靶标

功能	Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.
序列相似性	Belongs to the serpin family. Ov-serpin subfamily.
翻译后修饰	The signal sequence is not cleaved.
细胞定位	Cytoplasm. Secreted > extracellular space.

图片



Western blot - Anti-SerpinB2 antibody (ab105739)

Anti-SerpinB2 antibody (ab105739) at 1/500 dilution + HeLa whole cell lysates at 50 µg

Secondary

Goat-anti-Rabbit HRP at 1/5000 dilution

Developed using the ECL technique.

Predicted band size: 47 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors