

Product datasheet

Anti-Smallpox B5R antibody ab53309

1 图像

概述

<b>产品名称</b>	Anti-Smallpox B5R抗体
<b>描述</b>	兔多克隆抗体to Smallpox B5R
<b>经测试应用</b>	<b>适用于: WB</b>
<b>免疫原</b>	Synthetic peptide corresponding to Smallpox B5R aa 180-194. Variola virus Sequence: HYKKSESESYNELIK

 [Run BLAST with](#)

 [Run BLAST with](#)

性能

<b>形式</b>	Liquid
<b>存放说明</b>	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
<b>存储溶液</b>	Preservative: 0.05% Sodium azide Constituents: 0.2% Gelatin, PBS
<b>纯度</b>	Protein G purified
<b>克隆</b>	多克隆
<b>同种型</b>	IgG

应用

Our [Abpromise guarantee](#) covers the use of **ab53309** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 0.1 - 1 µg/ml. Detects a band of approximately 30 kDa (predicted molecular weight: 64 kDa).

靶标


**相关性** B5R is an EEV-specific Type I membrane glycoprotein essential for EEV formation. Vaccinia

virus, a member of Poxviridae family, is a dsDNA enveloped virus. It consists of 4 types of virion stages among which extracellular enveloped virus (EEV) is critical for cell-to-cell and long range viral infection. Some cellular antigens are associated with EEV and these vary with the cell type used to grow the virus. B5R consists of a large ectodomain containing regions similar to SCR domains of complement regulatory proteins followed by a transmembrane region and a cytoplasmic tail. B5R is essential for efficient wrapping of IMV, actin tail formation, normal plaque size, virus virulence, and for EEV formation.

细胞定位

Cell Membrane

### Anti-Smallpox B5R antibody 图像



**All lanes :** Anti-Smallpox B5R antibody (ab53309) at 0.1 µg/ml

**Lane 1 :** recombinant fusion protein containing amino acids 180-194

**Lane 2 :** fusion partner only

**Predicted band size :** 64 kDa

**Observed band size :** 30,60 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors