

## Product datasheet

# Anti-Staphylococcus DNase antibody ab15945

### 1 References

#### 概述

产品名称	Anti-Staphylococcus DNase抗体
描述	羊多克隆抗体to Staphylococcus DNase
宿主	Sheep
经测试应用	适用于: WB, ELISA
种属反应性	与反应: Staphylococcus aureus
免疫原	Full length native protein (purified).

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Upon delivery aliquot. Avoid freeze / thaw cycle.
存储溶液	Preservative: None Constituents: 0.9% Sodium chloride, 0.01M PBS, pH 7.4
纯度	Immunogen affinity purified
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab15945** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
WB		Use a concentration of 1 µg/ml.
ELISA		Use at an assay dependent dilution.

#### 靶标

## 相关性

One of the important virulence factors in *S. aureus* is DNase. This can be used to distinguish *Staphylococcus aureus* from the coagulase-negative staphylococci by the so called DNase test.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors