

## Product datasheet

# Anti-Staphylococcus DNase antibody ab15945

### 1 References

#### 概述

|       |  |
|-------|--|
| 产品名称  | Anti-Staphylococcus DNase抗体            |
| 描述    | 羊多克隆抗体to Staphylococcus DNase          |
| 经测试应用 | 适用于: WB, ELISA                         |
| 种属反应性 | 与反应: Staphylococcus aureus             |
| 免疫原   | Full length native protein (purified). |

#### 性能

|      |   |
|------|---|
| 形式   | Liquid  |
| 存放说明 | Shipped at 4°C. Upon delivery aliquot. Avoid freeze / thaw cycle.           |
| 存储溶液 | Preservative: None<br>Constituents: 0.9% Sodium chloride, 0.01M PBS, pH 7.4 |
| 纯度   | Immunogen affinity purified   |
| 克隆   | 多克隆   |
| 同种型  | IgG   |

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab15945** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| 应用    | Ab评论 | 说明                                  |
|-------|------|-------------------------------------|
| WB    |      | Use a concentration of 1 µg/ml.     |
| ELISA |      | Use at an assay dependent dilution. |

#### 靶标

**相关性** One of the important virulence factors in *S. aureus* is DNase. This can be used to distinguish *Staphylococcus aureus* from the coagulase-negative staphylococci by the so called DNase test.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### **Our Abpromise to you: Quality guaranteed and expert technical support**

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

### **Terms and conditions**

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors