

## Product datasheet

# Anti-Staphylococcus aureus antibody (FITC) ab68950

### 1 References

#### 概述

产品名称	Anti-Staphylococcus aureus抗体(FITC)
描述	兔多克隆抗体to Staphylococcus aureus (FITC)
偶联物	FITC. Ex: 493nm, Em: 528nm
经测试应用	适用于: ICC/IF, IHC-Fr
种属反应性	与反应: Staphylococcus aureus
免疫原	Soluble and structural antigens of Staphylococcus aureus whole bacterium; ATCC #27660.
常规说明	This antibody reacts with Staphylococcus aureus in bacterial and infected tissue samples.

#### 性能

形式	Liquid
存放说明	Shipped at 4°C. Store at +4°C.
存储溶液	Preservative: 0.1% Sodium Azide Constituents: PBS, 10mg/ml BSA
纯度	Whole antiserum
Primary antibody说明	This antibody reacts with Staphylococcus aureus in bacterial and infected tissue samples.
克隆	多克隆
同种型	IgG

#### 应用

Our [Abpromise guarantee](#) covers the use of **ab68950** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

应用	Ab评论	说明
ICC/IF		1/10.
IHC-Fr		1/10.

## 靶标

---

**相关性** Staphylococcus aureus are gram-positive bacteria that, although part of the normal flora of humans (nasal passages, skin and mucous membranes), have the potential to become opportunistic pathogens. S. aureus is a frequent cause of a wide range of infections and some bacterial strains have become resistant to many commonly used antibiotics.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
  
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.cn/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors