

Streptavidin HRP (ready-to-use) ab64269

★★★★★ [1 Abreviews](#) [42 References](#) [2 图像](#)

概述

产品名称	Streptavidin HRP (ready-to-use)
经测试应用	适用于: IHC-P
常规说明	Streptavidin-HRP for use with biotinylated secondary antibodies during IHC / immunohistochemistry. Optimized, ready-to-use format. This product was previously called Streptavidin-Peroxidase.

IHC protocol suitable for use with Streptavidin-HRP ab64269:

1. Deparaffinize and rehydrate formalin-fixed paraffin-embedded tissue section.
2. Use appropriate **antigen retrieval buffer or enzyme** (primary antibody dependent) to treat sections. Wash 3 times in buffer.
3. Add enough **hydrogen peroxide blocking solution** to cover the sections. Incubate for 10 minutes. Wash 2 times in buffer. If necessary, use **avidin biotin blocking**.
4. Apply **protein block** (or **normal serum** from same species as secondary antibody) and incubate for 5 minutes at room temperature to block nonspecific background staining. Wash once in buffer.
5. Apply primary antibody in **antibody diluent** and incubate.
6. Wash 4 times in buffer. Incubate slide with **biotinylated secondary antibody** (or **HRP polymer secondary antibody** and skip step 7). Wash 4 times in buffer.
7. **Apply streptavidin-HRP ab64269 and incubate for 10 minutes at room temperature.**
8. Rinse 4 times in buffer. Place slide in **DAB substrate** or **AEC substrate** and incubate until desired color is achieved. Rinse 4 times in buffer.
9. Add enough drops of **hematoxylin** to cover the section. Incubate for 1 minute.
10. Rinse 7-8 times in tap water. Add **mounting medium** to cover the section.

Find complete IHC kits, and reagents for antigen retrieval, blocking, signal amplification, visualization, counterstaining, and mounting in the [IHC kits and reagents guide](#).

性能

形式	Liquid
----	--------

存放说明 Store at +4°C.

存储溶液 pH: 7.60
Preservative: 0.2% Proclin 300
Constituent: PBS

相关性 Streptavidin peroxidase is used in a labeled streptavidin-biotin immunoenzymatic antigen detection system. This technique involves the sequential incubation of the specimen with an unconjugated primary antibody specific to the target antigen, a biotinylated secondary antibody which reacts with the primary antibody, enzyme-labeled streptavidin, and substrate-chromogen.


应用

The Abpromise guarantee **Abpromise™**承诺保证使用ab64269于以下的经测试应用

“应用说明”部分 下显示的仅为推荐的起始稀释度;实际最佳的稀释度/浓度应由使用者检定。


应用	Ab评论	说明
IHC-P		Use at an assay dependent concentration.

图片



ab64269 Straptavidin Peroxidase used along with an anti-CD20 antibody on Human Tonsil

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Streptavidin Peroxidase (Ready to Use) (ab64269)



Human colon carcinoma stained with anti-p53 antibody using Streptavidin Peroxidase (ab64269).

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Streptavidin Peroxidase (Ready to Use) (ab64269)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.cn/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors